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Product warranty and after-sales are valid in japan.

Sensitivity test (P.3)

To prevent skin trouble, carry out a patch test in advance.

- When you use a cosmetic product for the first time, perform a patch test (skin test) and check your skin.
 - Wash the inside of your arm with soap, and wipe it down with a towel, or the like.
 - Apply the cosmetic product to the skin area you have washed and leave it as is for 48 hours.
 - Apply a 10-yen-coin-size light coating, and let it dry naturally.
 - Do not run water over or rub the patch-tested part.

CAUTION

If any abnormality occurs on the application part, such as a rash, reddening, itching or irritation, immediately stop testing and wash away.

 - Check that there is no abnormality.

The skin lotion can be used if no abnormality occurs within 48 hours.
 - If you find something unusual with your skin after the patch test, do not use the cosmetic product with this device.
- Check your skin by testing if this device and your regular cosmetic product go well together.
 - Wash the inside of your arm with soap and wipe dry with a towel.
 - Place this device on the skin area you have washed for about 10 seconds by following the directions on how to use "Step 2 Hot Moisturize," on Page 17.
 - You can use the device if there is nothing unusual with your skin after following the procedures in the above-mentioned patch test and leaving the skin as it is for 48 hours.

Preparation before use (P.4-7)

How to charge (P.4)


- Since this is not fully initially charged, charge before first using.
- Make sure to use only 100-240V and use the included power adaptor.
- When using the device overseas, check the plug with the "Guidelines for use overseas" on page 27.

 WARNING	Do not insert and remove the power plug from an electric outlet with wet hands. (Cause of electric shocks and short-circuits)
	Do not use in places splashed with water, such as in a bath, or in damp places. (Cause of failures or short-circuits)

- Turn off the switch, and connect the power adaptor to the main body and electric outlet.

- After the standard charging time (1.5 hour) has elapsed, pull out the power adaptor from the electric outlet and the main body.

- This device is for rechargeable use only. It cannot be powered on with the power adaptor connected.
- The main body and the power adaptor become slightly hot, however, this is not unusual.

Charge indicator lamp
During charge : The lamp turns on. (Red)
After charging is completed :
The lamp turns off.

- Charging time is 1.5 hour.
- The device can be used approximately 7 times on a single charge.
(When all modes are used at "High power")
Depending on the operating conditions such as the ambient temperature, operating times may decrease.

Time to charge (P.5)

If the charge indicator lamp flashes when the power switch is "On," it indicates that a charge is required. Charge as soon as possible.

- In order to make full use of battery performance, observe the following.
- If you have not used this product for three months or longer, charge before use.
 - Do not charge every time you use this product.
 - The recommended charging temperature is in the range of 10 - 35°C. The battery may take longer to charge or may not charge at all at temperatures outside of the recommended temperature range.
 - Do not charge with a power supply other than that for home use.
 - If any TV noise occurs while the TV is on during charge or use of the product while plugged into the same outlet, change to another outlet.

Before use (P.6)

 CAUTION	Never use the device when you have a wound or an inflammation on your skin. (Cause of problems with or damage to your skin.)
	Remove jewelry or contact lenses. (Cause of problems with your skin or poor health)
	Check the heating dimple head or cooling head for any damage or deformation before use. (Cause of damage to skin)

1. About mode switching

Mode	Cotton and Cosmetics
Hot Cleanse (cleans fine residue remaining in pores)	Use new cotton saturated with your regular skin lotion.
Hot Moisturize (performs moisturizing care by allowing moisturizing ingredients to penetrate into skin)	Use new cotton saturated with your regular skin lotion. If you use a gel-like or thick skin lotion, you can also apply it to your skin without using cotton.
Hot Mask (enhances the effect of milky lotion or sheet mask)	Use your regular milky lotion or a sheet mask by applying to skin. Do not use cotton.
Cool (tightens the pores of your skin by cooling effect)	Apply the cooling head directly to your skin. Do not use cotton and skin lotion.

2 Clean skin by washing your face so that makeup and dirt will wash off.

Wipe off makeup with a cleansing product.

3 Attach cotton to the body of the device.

(1) Remove the protection seal from the heating dimple head.

Unless the protection seal is removed, the device will not operate.


(2) Place cotton over the head and attach it with the included ring.

Use commercially available cotton. However, small cotton or thick cotton may not work in some cases.

The ring has an orientation with a clear direction. Attach it from straight above with the wider part of the inner circumference placed toward the heating dimple head. When it doesn't attached easily, please attach it while pressing right and left side of the ring.

Do not attach the ring in the reverse direction. Please be careful not to squeeze it, otherwise, it cannot be removed.

The narrower side should be upward.

4 How to remove the ring.


● Remove the ring by pulling it straight up while holding down the lower part.

5 Please observe the following so as not to damage your skin when you use the device.

(1) The maximum amount of time allowed for skin care with use of this device is 8 minutes a day per person in basic mode.

● Limit use of step 1, 2, 3 to no more than 3 minutes per day.

● Do not use step 4 continuously.

● Limit use of the device to no more than 11 minutes per day.

● It is recommended that you use "Hot Mask" once a day for milky lotion or once or twice a week for sheet mask.

(2) Use the device on bare skin after removing makeup.

(3) Do not use the device on the same area of skin for more than 10 seconds.

(4) Switch between high and low by pressing the level switch button according to the condition of your skin.

● You can choose high, medium and low as you like. If your skin tingles, use the device on "Low."

● If you feel something unusual with your skin, stop using the device and consult a skin care specialist.

● The device is equipped with a memory function to memorize the levels of "High," "Medium," and "Low" that you would normally use. For details, please see "About memory function" on Page 26.

Names of parts and their functions (P.8)

① Heating dimple head (titanium material)

The device vibrates while being used in modes other than "Cool."

② Mode indicator lamp

Mode	During preparation	During use
Hot Cleanse	Orange lamp flashing	Orange lamp lights
Hot Moisturize		
Hot Mask		
Cool	Blue lamp flashing	Blue lamp lights

③ Mode switch button

Every time the button is pressed, the mode for use will be switched.

④ Power switch (Off/On)

Turn on the power switch, and it will turn "On." Press it again, and it will turn "Off."

⑤ Body

⑥ Body socket

⑦ Level indicator lamp

Low (green): Medium (green): High (red)
When the power switch is switched to "On," the lamp indicating "High," "Medium," or "Low" is turned on.

- ⑧ Level switch button
Every time the button is pressed, the electrical potential will change.
There will be no change in vibration strength.
- ⑨ Charge indicator lamp
During charge : The lamp turns on. (Red)
After charging is completed : The lamp turns off.
- ⑩ Buzzer (Switching sound/ready sound/warning sound/abnormality detection sound/finish sound/indicating timer sound)
Every time the switch or button is pressed, the switching sound will beep. When the preparation is completed in "Cool" mode, the ready sound will beep. The warning sound will beep to indicate the continuous use of "Cool" mode. In the event of a failure, an abnormality detection sound will beep. When each mode is finished, the finish sound will beep.
- ⑪ Cooling head (aluminum material)
Apply this head to your skin in "Cool" mode.
- ⑫ Back plate (titanium material)
Touch this part and apply the heating dimple head to your skin in modes other than "Cool," and the device will start to operate.
It will not vibrate if you remove your hand from this plate.
- ⑬ Batteries (built-in)
Lithium ion battery 1 pcs.
- ⑭ Power adaptor
Model: KH-55
- ⑮ Power plug
- ⑯ Instrument plug
- ⑰ Cord
- ⑱ Accessories
- ⑲ Rings (two) : Rings for attaching cotton.
- ⑳ Protective cap

Safety precautions (P.9-12)

- Precautions described here are intended for safe and correct use of this product and are to prevent injury and damage to you and others before they occur.
- Precautions are classified as follows.
Make sure to observe the precautions because they are related to safety and important details.

 DANGER	Indicates that improper handling may result in death or serious injury to a user, and the extent of hazard that its urgency is high.
 WARNING	Indicates the extent of hazard that improper handling is predicted to result in death or serious injury to a user.
 CAUTION	Indicates the extent of hazard or damage that improper handling is predicted to result in slight injury to a user or damage to property (house, household effects, etc.).

<Example of pictorial indication>

	The  symbol indicates "Prohibited" (Do not do this).
	"Disassembly prohibited"
	"Use prohibited in wet areas"
	"Prohibited from being made wet"
	"Prohibited from use with wet hands"
	The  symbol indicates "Compulsory" (Make sure to do this).
	"General action"
	"Pull out the power cord from outlet"

DANGER

-  Persons who use the following medical electric devices should not use this product.
(This product may cause malfunction of medical electric devices.)
- Self-contained medical electric devices such as a pacemaker
 - Wearable medical electric devices such as an electrocardiograph

WARNING

-  Never use this product by other methods which are not described in this instruction manual.
(Cause of accidents, trouble and failure)
- Use this product only for face and neck. (Cause of accidents or trouble)
- The following persons should avoid use of this product. (Cause of accidents or trouble)
- Persons with a heart disease
 - Persons with a febrile disease
 - Persons with an infectious disease
 - Persons who have a malignant tumor
 - Women experiencing difficulty with pregnancy until right after birth
 - Persons who are menstruating
 - Individuals with pre-existing skin disorders such as atopic dermatitis or contact dermatitis
 - Persons with an allergic predisposition
 - Individuals who worry about skin during use
 - Persons under dental treatment
 - Persons whose bile capillaries have expanded due to long time use of steroid hormones or liver function disease
 - Persons who are under the influence of alcohol or medication, or who are very tired
 - Individuals with perception disorders as a result of circulatory failure due to diabetes
 - During serious acute inflammatory disorders (such as fatigue, shivering, blood pressure fluctuations)
 - Persons with an acute (pain) disease
 - Persons with a blood pressure disorder
 - Persons with a tubercular disease
 - Persons with facial neuralgia
 - Persons with hemophilia

Do not use this product on the following regions. (Cause of skin trouble etc.)

- Region where orthopedic treatment was given - Eye ball and upper eye lid
- Wound site - Facial melanoderma - Lips
- Region of mucosal membrane such as inside of mouth
- Region irritated due to pimples or skin rash
- Regions which suffer from algesthesia or perception disorder
- Regions where metal, plastic or silicone etc., is embedded under skin
- Regions where there is itching, redness, or psychotic fleck due to physical stimulus etc.
- Regions where skin irritation is caused by cosmetics

Use only 100-240V AC for power supply. (Cause of fires and electric shocks)

Keep away from small children. Do not use on the skin of a child. (Cause of accidents or failure)

Never allow pins or dust to adhere to the power plug and plug for equipment.

(Cause of electric shock, short circuit and ignition)

Never use a power cord or a power plug if they are damaged or insertion to the outlet is loosened.

(If the power cord breaks, fires or electric shocks may occur.)

Never damage the power cord, unnecessarily bend, pull, twist, bind, place a heavy object on it, pinch, process, or wind around the power adaptor, etc.

(If the power cord breaks, fires or electric shocks may occur.)

Do not use this product near inflammable substances (gasoline, benzene, thinner, etc.).

(Cause of explosion or fire)

Use only the attached power adaptor. And, do not charge any other product with the attached power adaptor.

(Possibility of electric shock and ignition)

When feeling pain or disorder during use of this product, and sensing anything wrong with your skin after use, stop using and consult a dermatologist. (Cause of skin trouble etc.)

If any unusual sound or abnormality occurs on the main body, immediately stop using and contact the distributor for repair together with the power adaptor. (Cause of accidents and trouble)

Follow the description of the instruction manual for the operating time and operation frequency.

(Cause of accidents and trouble)

Use this device at an ambient temperature from 15 to 35°C. (If you use "Cool" mode at low temperature, the head may be overcooled, damaging the skin.)

Individuals who use for the first time should use while checking their skin condition.

(Cause of skin trouble or poor health)

Do not use this product beyond the specified operating time. (Cause of skin trouble)

If any abnormality occurs during charge or use, immediately pull out the power plug of the power adaptor from the outlet. (Cause of accident or fire)

Make sure to disconnect the power cord from the outlet for maintenance.

Never connect or disconnect the power plug to/from the outlet with wet hands.

(Cause of electric shocks and short-circuits)

Absolutely never clean the power adaptor with water.

And, never soak it in water or splash water on it. (Cause of failures and short-circuits)

Never use, charge or leave in a place where water is splashed such as the bathroom or in a place with high humidity. (Cause of failures and short-circuits)

Never modify the shaver. Only a qualified technician may make repairs.

(Cause of fires, electric shocks, injuries and failures)

CAUTION

Check the heating dimple head or cooling head for any damage or deformation before use.

(Cause of damage to skin)

Never use this product together with other beauty treatment devices at the same time.

(Cause of failures and troubles)

Never use this product if there is some injury or wound to the skin. (Cause of skin trouble or pain)

Do not leave this product in a place where the temperature becomes high such as the inside of a vehicle. (Cause of failure)

Do not drop this product, or give a strong impact. (Cause of injury or failure)

Do not wipe with alcohol, benzene or thinner, and do not spray insecticide.

(Cause of cracking, electric shock and ignition)

Use only the attached power adaptor. And, do not charge any other product with the attached power adaptor.

(Possibility of electric shock and ignition)

Remove accessories such as earrings and necklaces, and contact lenses for use.

(Cause of accidents or trouble)

If there is any severe skin disorder or irregularity on the skin surface, or due to downy hair, in rare instances you may feel strong irritation. In this case, avoid using this product on that region.

(Cause of accidents or trouble)

Do not stop the cotton on one part of the skin but always keep it in motion.

(Cause of skin trouble or bad physical condition)

Pull out the power plug and device plug without holding the cord.

(Cause of breakage, electric shock, short-circuit and ignition)

When there are people around the device, pay attention to any sudden movement during use.

(Cause of injury)

When carrying in a bag or the like, place into a separate pouch. (Cause of damage to the main body or failure)

Make sure to pull out the power plug from the outlet unless used.

(Cause of electric shock and fire due to deterioration)

Example for use (P.13-14)

Press the mode switch button and set the mode of your preferred treatment.

Quick tighten treatment (Required time: 2 minutes)

- Step 4 Cool

When you want to tighten your pores or cool down your skin before applying make up, you can just use "Cool." After 10-second preparation for cooling, it only takes one minute and 50 seconds to tighten the pores of your skin. Press the mode switch button to set to "Cool."

Quick moisturize treatment (Required time: 3 minutes)

- Step 2 Hot Moisturize

When you do not have much time in the morning or when you need to go out suddenly, you can just use "Hot Moisturize." It only takes three minutes to moisturize your skin thoroughly.

Press the mode switch button to set to "Hot Moisturize."

Tighten treatment (Required time: 5 minutes)

- Step 1 Hot Cleanse

- Step 4 Cool

When you want to remove fine residue remaining in pores and tighten the pores of your skin, use "Cool" after using "Hot Cleanse." This treatment tightens the pores of your skin.

After the "Hot Cleanse" mode is finished, press the mode switch button to set to "Cool."

Easy moisture course (Required time 6 minutes)

- Step 1 Hot Cleanse

- Step 2 Hot Moisturize

Before you go out or when you have more time, use "Hot Moisturize" after using "Hot Cleanse."

This treatment moisturizes the clean skin thoroughly.

Daily moisture course (Basic mode) (Required time: 8 minutes)

- Step 1 Hot Cleanse

- Step 2 Hot Moisturize

- Step 4 Cool

Standard daily treatment

After using "Hot Cleanse," moisturize the clean skin in "Hot Moisturize and tighten the pores of your skin by cooling effect in "Cool." After the "Hot Moisturize" mode is finished, press the mode switch button to set to "Cool."

Daily full course (Required time: 11 minutes)

- Step 1 Hot Cleanse

- Step 2 Hot Moisturize

- Step 3 Hot Mask (use of milky lotion)

- Step 4 Cool

Care for your skin using all modes.

Fine vibrations in "Hot Mask" mode and gentle switching of the electrode will condition the skin.

The cooling effect of "Cool" mode will tighten the pores of your skin.

Weekly full course (Required time: 11 minutes)

- Step 1 Hot Cleanse

- Step 2 Hot Moisturize

- Step 3 Hot Mask (use of sheet mask)

- Step 4 Cool

Care for your skin once or twice a week by using all modes.

Fine vibrations in "Hot Mask" mode and gentle switching of the electrode will condition the skin.

The cooling effect of "Cool" mode will tighten the pores of your skin.

How to use (P.15-26)

This device will not operate with the power adaptor plugged into an outlet.

Remove the power plug of the power adaptor from the outlet to use.

- Use after removing the protective cap.
- The heating dimple head slowly warms after power is turned on.

 WARNING	Do not use the product around eyes and ears. (Cause of accidents or bad physical condition)
	If you feel something is unusual with your skin while using the device, stop using it immediately and consult a skin care specialist. (Cause of accident or poor health)

Step 1 Hot Cleanse (P.15)

This mode removes fine residue remaining in pores which cannot be removed by just washing.

1. Saturate cotton attached to the body with skin lotion after washing your face.
Please use your regular skin lotion or rubbing-off skin lotion with an exfoliating effect.
(Do not use milky lotion or a gel-like.)

- Saturate the cotton with lotion to the extent that lotion permeates the back of the cotton.

Some skin lotion (just with mineral spring or water) may be difficult to bring out a good effect.

2. Press the power switch.
The level indicator lamp will be turned on, and "Hot Cleanse" on the mode indicator lamp will start flashing. If a lamp in another mode is flashing, press the mode switch button to switch to "Hot Cleanse."
3. Choose "level" of electrical potential with the level switch button. Switch the "level" with the level switch button according to the condition of your skin.

- Every time the level indicator button is pressed, the high, medium, and low of the electrical potential will be switched. There will be no change in vibration strength.

It is recommended that when using the device for the first time or if you feel some discomfort, use the device in "Low."

4. Hold the device with your hand touching the back plate of the body.
 - Make sure to hold the body with your hand touching the back plate.
5. Apply cotton to the skin and move the device slowly.
 - The mode indicator lamp will switch and turn on, and the device will start to vibrate.

- The time for use is within 3 minutes.

Cautions

Unless you apply the cotton to your skin while touching the back plate, the device will not vibrate.

(skin sensing function)

If the heating dimple head is not vibrating, the device is not working.

If you leave the device as it is for 3 minutes without vibration, the device will automatically switch "Off."

(Automatic shutoff function)

Indicating timer function


Skin care time is indicated by the buzzer beeping about every 45 seconds as a guide.

[Example of skin care sequence]

Divide the area from your face to décolletage into four zones as shown in the illustration.

By moving the heating dimple head in the order of ①②Both cheeks→③T-zone→④décolletage, the whole area can be efficiently and uniformly cared for.

- Apply the cotton gently to the face and slowly move the device by sliding it from the center of the face to the outer side.

⚠ CAUTION

Do not keep the cotton on the same area of skin, but move it around.
(Cause of problems with your skin or poor health.)

- Add more skin lotion when there is no moisture in the cotton or when the device is difficult to slide.
 - Do not apply the device on the same area of skin for more than 10 seconds.
 - If you feel your skin tingling by micro current while using the device, press the level switch button and switch to "Low." If you still feel your skin tingling, use the device by adding more skin lotion.
6. After 3 minutes have elapsed, the device will automatically shut off. (Auto shutoff function)
 - The finish sound (beep beep) will inform you.
 7. Remove the ring and the cotton.
 - Wipe off the surface of the heating dimple head with a tissue.
 - Do not use the cotton again once it is used. It might cause problems with your skin.

Step 2 Hot Moisturize (P.17)


When using cotton

This mode cares for skin by allowing moisturizing ingredients of skin lotion to permeate the skin.

1. Attach new cotton to the body with the attachment ring.
 - For attaching cotton, see Page 4.
2. Saturate the cotton with your regular skin lotion.
 - Do not use milky lotion or a gel-like or thick skin lotion.
 - Saturate the cotton with lotion to the extent that lotion permeates the back of the cotton.

- 3 Press the power switch.
- "Low" on the level indicator lamp lights and "Hot Moisturize" on the mode indicator lamp starts flashing.
- If a lamp in another mode is flashing, press the mode switch button to switch to "Hot Moisturize."


When the "Hot Cleanse" mode is finished, press the power switch and the "Hot Moisturize" mode will automatically be set.

- 4 Choose "level" of electrical potential with the level switch button. Switch the "level" with the level switch button according to the condition of your skin.
- Every time the level indicator button is pressed, the high, medium, and low of the electrical potential will be switched. There will be no change in vibration strength.
 - The memory function memorizes the "level" when the power has been shut off.

It is recommended that when using the device for the first time or if you feel some discomfort, use the device in "Low."

- 5 Hold the device with your hand touching the back plate of the body.
- Make sure to hold the body with your hand touching the back plate.
- 6 Apply cotton to the skin and move the device slowly.
- The mode indicator lamp will switch and turn on, and the device will start to vibrate.
 - The time for use is within 3 minutes.

Cautions

Unless you apply the cotton to your skin while touching the back plate, the device will not vibrate.

(skin sensing function)

If the heating dimple head is not vibrating, the device is not working.

If you leave the device as it is for 3 minutes without vibration, the device will automatically switch "Off."

(Automatic shutoff function)

Indicating timer function

Skin care time is indicated by the buzzer beeping about every 45 seconds as a guide.

[Example of skin care sequence]

Divide the area from your face to décolletage into four zones as shown in the illustration.

By moving the heating dimple head in the order of ①②Both cheeks→③T-zone→④décolletage, the whole area can be efficiently and uniformly cared for.

- Apply the cotton gently to the face and slowly move the device by sliding it from the center of the face to the outer side.

CAUTION

Do not keep the cotton on the same area of skin, but move it around.
(Cause of problems with your skin or poor health.)

- Add more skin lotion when there is no moisture in the cotton or when the device is difficult to slide.
 - Do not apply the device on the same area of skin for more than 10 seconds.
 - If you feel your skin tingling by micro current while using the device, press the Level switch button and switch to "Low."
If you still feel your skin tingling, use the device by adding more skin lotion.
- 7 After 3 minutes have elapsed, the device will automatically shut off. (Auto shutoff function)
- The finish sound (beep beep) will inform you.
- 8 Remove the ring and the cotton.
- Wipe off the surface of the head with a tissue.
 - Do not use the cotton again once it is used. It might cause problems with your skin.


When not using cotton

If you use a gel-like cosmetic (all-in-one type skin lotion for moisturizing) which is difficult to permeate the cotton or thick skin lotion in "Moisturize," you can use this device without attaching cotton to the body.

- 1 Apply a gel-like cosmetic or thick skin lotion that you would normally use onto the heating dimple head and smoothen it over the entire head with your fingers.

- 2 After confirming that the power adaptor is connected to the main body, press the power switch.
- "Low" on the level indicator lamp lights and "Hot Moisturize" on the mode indicator lamp starts flashing. If a lamp in another mode is flashing, press the mode switch button to switch to "Hot Moisturize."

When the "Hot Cleanse" mode is finished, press the power switch and the "Hot Moisturize" mode will automatically be set.


- 3 Choose "level" of electrical potential with the level switch button. Switch the "level" with the Level switch button according to the condition of your skin.
- Every time the level indicator button is pressed, the high, medium, and low of the electrical potential will be switched. There will be no change in vibration strength.
 - The memory function memorizes the "level" when the power has been shut off.

It is recommended that when using the device for the first time or if you feel some discomfort, use the device in "Low."

- 4 Hold the device with your hand touching the back plate of the body.
- Make sure to hold the body with your hand touching the back plate.

- 5 Apply the heating dimple head directly to the skin and move it slowly.
- The mode indicator lamp will switch and turn on, and the device will start to vibrate.
 - The time for use is within 3 minutes.

Cautions

Unless you apply the cotton to your skin while touching the back plate, the device will not vibrate.

(skin sensing function)

If the heating dimple head is not vibrating, the device is not working.

If you leave the device as it is for 3 minutes without vibration, the device will automatically switch "Off."

(Automatic shutoff function)

Indicating timer function

Skin care time is indicated by the buzzer beeping about every 45 seconds as a guide.

[Example of skin care sequence]

Divide the area from your face to décolletage into four zones as shown in the illustration.

By moving the heating dimple head in the order of ①②Both cheeks→③T-zone→④décolletage, the whole area can be efficiently and uniformly cared for.

- Add more skin lotion or other cosmetics when the heating dimple head does not easily slide.

- Apply the heating dimple head gently to the face and slowly move the device by sliding it from the center of the face to the outer side.

⚠ CAUTION

Do not keep the heating dimple head in the same area of the skin, but move it around. (Cause of problems with your skin or poor health)

- Do not apply the device on the same area of skin for more than 10 seconds.
 - If you feel your skin tingling by micro current while using the device, press the Level switch button and switch to "Low." If you still feel your skin tingling, use the device by adding more skin lotion or other cosmetics.
- 6 After 3 minutes have elapsed, the device will automatically shut off. (Auto shutoff function)
- The finish sound (beep beep) will inform you.
- 7 Wipe off the cosmetic on the heating dimple head and the body of the device with a tissue.

Step 3 Hot Mask (P.22)

(It is recommended that you use "Mask" once a day for milky lotion or once or twice a week for sheet mask.)


Fine vibrations and gentle switching of the electrode enhance the effect of milky lotion or sheet mask.

- Place a sheet mask on your face or spread milky lotion or cream that you would normally use over your entire face.
 - Wipe off any cosmetic on your hands.
- Press the power switch.
 - "Low" on the level indicator lamp lights and "Hot Mask" on the mode indicator lamp starts flashing. If a lamp in another mode is flashing, press the mode switch button to switch to "Hot Mask."

2 Press the power switch.

"Low" on the level indicator lamp lights and "Hot Mask" on the mode indicator lamp starts flashing. If a lamp in another mode is flashing, press the mode switch button to switch to "Hot Mask."

When the "Hot Moisturize" mode is finished, press the power switch and the "Hot Mask" mode will automatically be set.

3 Choose "level" of electrical potential with the level switch button.

Switch the "level" with the level switch button according to the condition of your skin.

- Every time the level indicator button is pressed, the high, medium, and low of the electrical potential will be switched. There will be no change in vibration strength.

It is recommended that when using the device for the first time or if you feel some discomfort, use the device in "Low."

4 Hold the device with your hand touching the back plate of the body.

- Make sure to hold the body with your hand touching the back plate.

5 Apply and move the heating dimple head slowly.

When using a sheet mask, slide the head over the sheet mask.

- The mode indicator lamp will switch and turn on, and the device will start to vibrate.
- The time for use is within 3 minutes.

Cautions

Unless you apply the heating dimple head to your skin while touching the back plate, the device will not vibrate. (skin sensing function)

If the heating dimple head is not vibrating, the device is not working.

If you leave the device as it is for 3 minutes without vibration, the device will automatically switch "Off." (Automatic shutoff function)

Indicating timer function


Skin care time is indicated by the buzzer beeping about every 45 seconds as a guide.

[Example of skin care sequence]

Divide the area from your face to décolletage into four zones as shown in the illustration.

By moving the heating dimple head in the order of ①②③④ Both cheeks → ③ T-zone → ④ décolletage, the whole area can be efficiently and uniformly cared for.

- Apply the heating dimple head gently to the face and slowly move the device by sliding it from the center of the face to the outer side.

CAUTION

Do not keep the heating dimple head in the same area of the skin, but move it around. (Cause of problems with your skin or poor health)

- Do not apply the device on the same area of skin for more than 10 seconds.
- If you feel your skin tingling by micro current while using the device, press the Level switch button and switch to "Low."

6 After 5 minutes have elapsed, the device will automatically shut off. (Auto shutoff function)

- The finish sound (beep beep) will inform you.

7 Wipe off the cosmetic on the heating dimple head and the body of the device with a tissue.

Step 4 Cool (P.24)

In this mode, the head surface that is cooled about 10°C below the ambient temperature tightens the pores of your skin. The cool touch to the face is refreshing.

WARNING

Use this device at an ambient temperature from 15 to 35°C. (If you use "Cool" mode at low temperature, the head may be overcooled, damaging the skin.)

Cautions

- Continuous use of the "COOL" mode makes it difficult to cool the cooling head.


The warning sound and the flashing "COOL" on the mode indicator lamp will inform you. The warning sound can be canceled by pressing the power switch, but the "COOL" mode cannot be used for about one hour after cancellation.

"HOT CLEANSING," "HOT MOISTURIZE" and "HOT MASK" can still be used.

1 After confirming that the power adaptor is connected to the main body, press the power switch.

- The level indicator lamp will go out, and "Cool" on the mode indicator lamp starts flashing. If a lamp in another mode is flashing, press the mode switch button to switch to "Cool."

When the "Hot Mask" mode is finished, press the power switch and the "Cool" mode will automatically be set.

- 2 Wait until the cooling head has cooled.
- Preparation for cooling takes about 10 seconds. "Cool" on the mode indicator lamp flashes during preparation.
 - When the device is ready for use, the ready sound will beep, and the mode indicator lamp will be turned on.

Caution

The device does not vibrate in "Cool" mode.

- 3 Apply the cooling head gently to your skin and move it slowly just like patting your face.

- The mode indicator lamp will change to lighting from flashing.
- The cooling head does not vibrate.
- The time for use is within one minute and 50 seconds.
- Since you can apply the device to the same area of your skin for only one second, repeatedly apply the head to your skin on and off towards the edges of the face.
- There is no indicating timer function in the "COOL" mode.

Caution

- Do not slide the cooling head over the skin.
 - Do not keep the cooling head in the same area of the skin, but move it around.
- The cooling effect of the cooling head may be lessened.

- 4 After one minute and 50 seconds have elapsed, the device will be automatically shut off. (Automatic shutoff function)

- The finish sound (beep beep) will inform you.
- When the "Cool" mode is finished, the main body and the cooling head become slightly hot, however this is not unusual.

Caution

If all of the mode indicator lamp and the level indicator lamp start flashing and the abnormality detection sound beeps during use, this is a warning of device failure.

Pull out the power adaptor from the electric outlet, and contact the distributor for repair.

About skin sensing function (P.26)

- This device is a product which operates through the human body.
- Applying the skin to the side plate and the heating dimple head will make an electrical circuit start to operate. Unless you touch the side plate and apply the heating dimple head to the skin, it will not start to operate.

About memory function (P.26)

- This device is equipped with a memory function to store memory of the electrical potential levels of "High," "Medium," or "Low".

In each usage mode (Step 1 to 3), the device stores memory of the level just before the automatic shut off or the power switch is turned "off."

If the level switch button is pressed while the device is being used, the electric potential "level" can be switched.

Care for main body (P.28)

- Wipe away makeup attached to the heating dimple head, cooling head, main body and ring with tissue paper or a cloth.
- When dirt is difficult to remove, wipe it off with a cloth after saturating it with water and wring it thoroughly.

⚠ WARNING	Do not clean the main body with water. Or do not dip in water, and do not pour water. (Cause of failure)
⚠ CAUTION	Make sure to turn OFF the switch. (Improperly turning on the switch may give a tingling sensation or result in malfunction or short-circuit.)
	Do not use alcohol, benzene or thinner. (Cause of deformation or discoloration)

How to preserve (P.28)

- Attach the included protective cap to the main body and store it in a place with little humidity and dust after cleaning it.
- The protective cap can be attached with the ring attached to the main body.

⚠ WARNING	Do not leave in a water splashing area or an area with significant humidity, such as a bath. (Cause of short-circuits or failure)
	Do not wind the power cord around the power adaptor or the main body. (Cause of damage to cord)

Cautions

Do not leave in an area subject to direct sunlight (inside of a vehicle in summer, etc.) or near heaters. (Cause of deformation or discoloration)

- Align the hook with the neck to attach.

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